

# EasyOne Air

All the Portable Advantages,  
One Connected Solution

NicheMedical  
Solutions in Respiratory Care & Diagnosis



## Spirometry (FVC, FVL, Tidal FVC, Tidal FVL, SVC & MVV)

The proven ultrasound technology  
**NDD TrueFlow**

**no calibration, no warm-up  
time, no moving parts**

Automated user guidance throughout maneuvers based on  
ATS/ERS standards 2019 and 2005

Quickly assess test quality with full color real time  
curves and instant interpretation

Data exchange via Bluetooth to the PC

Intuitive user guidance

Large color touch screen for easy data entry and navigation

Pediatric incentive via Bluetooth to EasyOne Connect

EasyOne Connect for seamless EMR integration

Rechargeable battery

**NDD TrueFlow**  
makes the difference

The original ultrasonic flow measurement is highly accurate in all flow ranges, independent of gas composition, pressure, temperature and humidity and does not require calibration during its lifetime. The sensor is never in direct contact with the patient's flow. NDD TrueFlow is a hygienic and resistance-free solution.

**NDD EasyConnect**  
intelligent interfacing

NDD's connectivity engine offers a comprehensive set of default configured HL7 and XML interfaces. With one database and one platform for all EasyOne point-of-care solutions, data management has never been easier.

### Standards & Recommendations

#### Quality, Medical Devices & Electrical

ISO 13485, ISO 14971, IEC 62366,  
IEC 62304, ISO 26782, ISO 23747,  
IEC 60601-1, IEC 60601-1-2, ISO  
10993-1

#### FDA

510(k) clearance

#### Associations & Institutes

ATS/ERS 2019 and 2005, NIOSH, OSHA,  
SSA Disability

### Languages

English, Danish, Dutch, French, German, Italian, Polish, Portuguese,  
Russian, Spanish, Swedish

### Technical

#### Printing options

Direct to printer or with EasyOne Connect  
software

#### Data management

EasyOne Connect

#### Export

HL7, XML, GDT, with software

#### Data links

USB, Bluetooth

#### Test storage

Up to 10'000 tests

#### Age range

Spirometry > 4 years

#### Dimensions

87 x 155 x 36 mm (H x B x T), 356 g  
3.4 x 6.1 x 1.4" (H x W x D), 13 oz

#### Device classification

Type BF applied part

#### Operating conditions

Temp 0 - 40 °C / 32 - 104 °F  
Rel. Humidity 5 - 90%  
Atmosph. Pressure 700 - 1060 hPa

#### Power supply

5 VDC, Standby 0.3W

#### Rechargeable battery

Exchangeable, 3.6 VDC

## Parameters

<b>FVC</b>	ATI, BEV, EOTV, FEF10, FEF25, FEF25-75, FEF25-75_6, FEF40, FEF50, FEF50/FVC, FEF50/VCmax, FEF60, FEF75, FEF75-85, FEF80, FET, FET25-75, FEV.25, FEV.5, FEV.5/FVC, FEV.75, FEV.75/FEV6, FEV.75/FVC, FEV.75/VCmax, FEV1, FEV1/FEV6, FEV1/FVC, FEV1/FVC6, FEV1/VC, FEV1/VCmax, FEV3/FVC, FEV3/VCmax, FEV3, FEV6, FVC, MEF20, MEF25, MEF40, MEF50, MEF60, MEF75, MEF90, MMEF, MTC1, MTC2, MTC3, MTCR, PEF, PEFT, t0, VC, VCmax
<b>FVL</b>	ATI, BEV, CVI, E50/150, EOTV, FEF10, FEF25, FEF25-75, FEF25-75_6, FEF40, FEF50, FEF50/FVC, FEF50/VCmax, FEF60, FEF75, FEF75-85, FEF80, FET, FET25-75, FEV.25, FEV.5, FEV.5/FVC, FEV.75, FEV.75/FEV6, FEV.75/FVC, FEV.75/VCmax, FEV1, FEV1/FEV6, FEV1/FIV1, FEV1/FIVC, FEV1/FVC, FEV1/VC, FEV1/VCmax, FEV3/FVC, FEV3/VCmax, FEV3, FEV6, FIF25, FIF25-75, FIF50, FIF50/FEF50, FIF75, FIV.25, FIV.5, FIV1, FIVC, FVC, MEF20, MEF25, MEF40, MEF50, MEF60, MEF75, MEF90, MIF25, MIF50, MIF75, MMEF, MMIF, MTC1, MTC2, MTC3, MTCR, PEF, PEFT, PIF, t0, VC, VCmax
<b>SVC</b>	ERV, IC, IRV, Rf, VC, VCex, VCin, VCmax, VT
<b>MVV</b>	MVV, MVV6, MVVtime, Rf, VCext, VT

## Predicted normal values Spirometry

<b>GLI</b>	Quanjer 2012, Stanojevic 2009
<b>North America</b>	Crapo 1981, Dockery (Harvard) 1993, Eigen 2001, Gutierrez (Canada) 2004, Hsu 1979, Knudson 1983, Knudson 1976, Morris 1971 & 1976, NHANES III (Hankinson) 1999, Polgar 1971
<b>Latin America</b>	Chile 2010, Chile (Pediatrics) 1997, Pereira 1992, Pereira 2006/2008, Pérez-Padilla (PLATINO) 2006, Pérez-Padilla (Mexico) 2001, Pérez-Padilla (Mexico, Pediatrics) 2003
<b>Europe</b>	ERS (ECCS, EGKS, Quanjer) 1993, Garcia-Rio (SEPAR) 2013, Falaschetti 2004, Forche (Austria) 1988 & 1994, Klement (Russia) 1986, Roca (Spain, SEPAR) 1982, Rosenthal 1993, Sapaldia (Switzerland) 1996, Vilozni 2005, Zapletal 1977, Zapletal 2003
<b>Europe Scandinavia</b>	Berglund Birath (Sweden) 1963, Finnish 1982 (1998), Gulsvik (Norway) 1985, Hedenström 1985 & 1986, Langhammer (Norway) 2001, Kainu (Finland), 2016, Nystad 2002
<b>Australia</b>	Gore Crockett 1995, Hibbert 1989
<b>Asia</b>	Chhabra (India) 2014, Dejsomritrutai (Thailand) 2000, Indonesia 1992, IP (China, HongKong) 2000 & 2006, JRS 2001 & 2014
<b>Africa</b>	Mengesha (Ethiopia), 1985

## Flow/Volume Sensor

<b>Type</b>	Ultrasonic transit time
<b>Flow Range</b>	± 16 l/s
<b>Flow Resolution</b>	4 ml/s
<b>Flow Accuracy (except PEF)</b>	± 2% or 0.020 l/s
<b>PEF Accuracy</b>	± 5% or 0.200 l/s
<b>MVV Accuracy</b>	± 5% or 5 l/min
<b>Volume Range</b>	± 12 l
<b>Volume Resolution</b>	1 ml
<b>Volume Accuracy</b>	± 2% or 0.050 l
<b>Resistance</b>	0.3 cm H <sub>2</sub> O/l/s at 16 l/s

**EasyOne Air** complete package with device, power plug and adapter, cradle, USB cables, battery pack etc.

## Order Information

Part Number	Product
2500-2A	EasyOne Air US
2500-2INT	EasyOne Air International

## Accessories

Part Number	Product
5050-50	EasyOne FlowTube, individually wrapped, box of 50 pcs
5050-200	EasyOne FlowTube, individually wrapped, box of 200 pcs
5050-500	EasyOne FlowTube, individually wrapped, box of 500 pcs
2030-2	NDD Calibration syringe 3L with EasyOne FlowTube Cal Check Adapter
5030-2	EasyOne FlowTube Adapter
2500-50.1	EasyOne Air USB cable B-micro (cradle to printer)
2500-50.2	EasyOne Air bluetooth dongle
2500-50.4	EasyOne Air battery pack
2500-50.5	EasyOne Air power supply with adapters
2500-50.11	EasyOne Air cradle with cables, power supply and adapters