

SideStream Plus

Breath enhanced nebuliser

SideStream Plus – reusable nebuliser, single patient use

Product description

	Manufacturer part number	Unit of issue	Replacement for VentStream part number	NHS Supply Chain code	SideStream Plus handset	Valved mouthpiece	Air supply tubing	Filter holder plus 10 filters	Filters (x 50)	Valved tee piece	Non-valved mouthpiece	Wide bore aerosol hose
SS+ with expiratory filter holder & filters	1091999	Eaches	1170A 1170	FDE039	•	•	•	•		•	•	
SS+ with aerosol ventilation hose (1m)	1092003	Eaches	1171A 1171	FDE038	•	•	•			•	•	•
SS+ with valved tee piece & mouthpieces	1092118	Eaches	1182A 1183	–	•	•	•			•	•	
Expiratory filters	1093190	50 pack	1220	FTC155					•			
Expiratory filter holder & 10 filters	1091917	Eaches	1221	FTC210				•		•	•	
Aerosol ventilation hose (1m)	1093185	Eaches	1080A	–								•



SS+ with expiratory filter holder & filters



SS+ with aerosol ventilation hose



SS+ with valved tee piece & mouthpieces

- If a mask is required for nebuliser treatment:
 - replace the mouthpiece with a sealing mask
 - ensure a valved tee piece is connected to the sealing mask
- SideStream Plus requires a compressor driving flow of 3.5 LPM minimum (suitable for use with Porta-neb and InspirationElite)
- Refer to the 'Instructions For Use' for all product information