

spirette™

Breathing Tube for
Infection Control



**Avoid cross-contamination
and reduce cleaning to an
absolute minimum with the
single-patient-use spirette.**



All parts exposed to patient's breath
are **single-use** and **disposable**

No expiration date

No special storage conditions

Technical data

Material

spirette: Linear low density polyethylen (LLDPE)
Bag: Polyethylen (PE)

Transport and storage information

Temp. -20 to 50°C
Rel. Humidity 5 to 95%
Athmosph. Pressure 400 to 1060 hPa
no expiration date

MDD 93/42/EEC

Class I

Dimension (single spirette)

25 x 38 x 140 mm (H x W x D), 8 g

Compatible with

EasyOne, Easy on-PC, EasyOne Pro and EasyOne Pro LAB

Order information

Part Number

Description

2050-1

spirette, individually wrapped, dispenser box of 50 pcs

2050-5

spirette, individually wrapped, box of 200 pcs (dispenser box included)

2050-10

spirette, individually wrapped, box of 500 pcs (without dispenser box)

NicheMedical
Solutions in Respiratory Care & Diagnosis

Tel: 1300 136 855 Email: info@nichemedical.com.au
www.nichemedical.com.au